Advantages and benefits

- Ultimativ SAXS/WAXS/USAXS/GISAXS and GIWAXS Beamline for the lab
- self-developed sample environments integrable
- Variable size of incident X-ray beam adapted to requirements
- Screening with small beam diameter to analyze differences in a sample
- Statistical information about the size distribution of the particles
- Measurement of samples in the desired sample environment
- Statistically reliable information in the range of a few nm compared to imaging methods

Equipment

X-ray source	High-power micro-focus rotating-anode X-ray source Brilliant micro-focus X-ray source (air cooled)
Optics	Parabolic or hyperbolic multilayer
Collimation	Scatterless and automated slits
Linear motion	3 linear axes each for detector and sample positioning, additionally a 360° rotation axis for sample rotatability
Sample environment	Sample holder for up to 8 solid samples, capillary holder for up to 10 capillaries Optional holder with heater for 1 capillary (temperature range: 5 – 80 °C) Optional holder for powders with defined thickness
Detector	SAXS/WAXS Detector Dectris EIGER 1M (1035 x 1065 Pixel, Pixelsize: 75 µm) Optional WAXS Detector

Specifications

Sample detector Distance	Stepless between 6 cm and 3.2 m
Used wavelength	1.5 Å (Cu-Kα1)
Q-Scale	0.005 – 5 Å ⁻¹
Associated quantities in real space	5 Å – 150 nm
Optional USAXS Beamsize	Down to 2 x 10 ⁻⁵ Å ⁻¹ Between 300 µm and 1 mm (smaller diameters possible with other optics)

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Kontakt
Dr. Bernhard Schummer
bernhard.schummer@iis.fraunhofer.de

Nanostrukturanalyse/
Nanostructural analysis

Flexible SAXS/WAXS/USAXS/GISAXS
and GIWAXS Beamline

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

Plot of the reduced scattering profiles of plant starch dissolved in water. Each of the curves shown is composed of four individual measurements. The total measurement time for one curve is 8 hours.

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Size distribution of nanodiamonds separated by different sorting methods (Treatment 1-3) to obtain small particles.

Motivation

Small-angle X-ray scattering (SAXS) is a common method for material characterization. Possible samples can be investigated both in solid phase (e.g. powders, alloys, aerogels, wood, collagen) and as solution (e.g. colloidal solutions, polymer solutions, stabilized nanoparticle solutions, drug formulations). This scattering method finds its application not only in basic research and development, but also in quality control during in-line testing.

Custom Setup

With the current SAXS/WAXS setup, the distance between sample and detector can be continuously adjusted between 6 cm and 3.5 m. This allows a combination of small and wide angle measurements without sample change and changing the setup. The resulting reduced scattering profiles provide extremely high scattering intensities and over three full orders of magnitude ($0.005 - 5 \text{ \AA}^{-1}$) in momentum transfer Q for suitable samples (Fig. 1). With a USAXS module the measuring range can be extended down to $2 \times 10^{-5} \text{ \AA}^{-1}$. Without modification the sample can be measured


In principle, SAXS can be used to obtain information about the

- size and shape of the components of the sample.
- structure of the sample.
- orientation of the particles.

In contrast to electron microscopy, the results obtained in this way are not only valid for a small sample. By averaging over a very large number of particles in the sample, statistically relevant information is possible, e.g. for the size distribution.

in GISAXS or GIWAXS geometry.

The sample chamber is designed to be flexible enough to accommodate a wide variety of sample environments. Since the software is open, custom sample environments developed by the customer can also be integrated and automated. In addition, the measured intensity can be calibrated to match the differential cross section of the underlying scattering process. Such calibration can be used to determine more complex properties, such as the electron density distribution in the sample or the hydration level of colloids in solution.

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Scattering signal for a 30 weight percent (wtp) aqueous Pluronic P123 solution without a) and with b) stabilized cadmium sulfide nanoparticles at a temperature of 30°C (exposure time: two hours in each case). In a), the liquid crystalline phase of the micelles can be clearly seen due to the Bragg reflections that appear. The stabilized cadmium sulfide nanoparticles in b) drastically changes the detector image.

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Titelbild: ©Universität Würzburg